

Coronis 2MP (MFGD-2320)

2 MegaPixel premium diagnostic display system



PRODUCT SPECIFICATIONS**CORONIS 2MP (MFGD-2320)**

Analog	VGA D-SUB15 connector, BNC
Approvals	Safety: <ul style="list-style-type: none">■ CE, UL60950, CSA C22.2 No 60950 (c-UL), CB, IEC 60950, DEMKO EN 60950 Emission: <ul style="list-style-type: none">■ EN 55022, EN 50082-1, FDA 510 k
Color imaging	No
Dark room contrast	1000:1 typical
Display technology	<ul style="list-style-type: none">■ TFT AMLCD VA Mode LCD■ Grayscale, 3 subpixels per pixel
DVI	Complying to VESA DVI
Ergonomics	<ul style="list-style-type: none">■ Tilt & swivel base with P/L rotation■ 100 mm VESA mounting standard■ AR coated front cover
Gray imaging	True
Luminance	400 cd/m ² (116.7 fL) calibrated, 600 cd/m ² (175.1 fL)
Resolution	2MP (1600 x 1200 / 1200 x 1600)
Stabilization, calibration and automated QA	Embedded I-Guard
Supported resolution	<ul style="list-style-type: none">■ 1600 x 1200 landscape■ 1200 x 1600 portrait (Optional portrait Accelerator)■ VGA boot
USB	USB-hub with 1 up-and 2 downstream ports
Viewing angle	170° (at 10:1 contrast)

Last updated: 22 Oct 2018

Technical specifications are subject to change without prior notice. Please check www.barco.com for the latest information.